

# NETZSCH

## Adiabatic Reaction Calorimetry



Leading Thermal Analysis. ■

MMC 274<sup>™</sup>  
*Nexus*

# Multi Modular Adiabatic Calorimeter

## MMC 274 *Nexus*<sup>TM</sup>

The MMC 274 *Nexus*<sup>TM</sup> can be used in conjunction with the DSC and the Accelerating Rate Calorimeter where its special features better match the application. The MMC 274 *Nexus*<sup>TM</sup> is used for measuring heats of reaction, reaction rates, endotherms, heat capacities, phase changes, gas generation rates and vapor pressures.

### Like a DSC, the MMC 274 *Nexus*<sup>TM</sup> can:

- Measure sample power output as a function of temperature
- Measure endothermic and exothermic events
- Run in scanning modes or isothermal mode
- Utilize short test time
- Measure heat capacities

### Like an Accelerating Rate Calorimeter, the MMC 274 *Nexus*<sup>TM</sup> allows:

- Mixing, injection, venting
- Collection of pressure data as a function of sample temperature
- Adiabatic operation, including low thermal inertia tests
- Working with larger multi-phase samples



### Modular System

The MMC 274 *Nexus*<sup>TM</sup> consists of two parts: the base unit with electronics and the exchangeable calorimeter module. This guarantees maximum flexibility.

Exchangeable calorimeter modules are available for different applications:

- Accelerating Rate Calorimeter type module for safety testing
- *VariPhi*<sup>TM</sup> module for fire exposure and other tests
- External sample heater module for scanning tests
- Special modules for battery testing

### Multiple Testing Modes

The MMC 274 *Nexus*<sup>TM</sup> can work in different testing modes. The scanning methods provide sensitivities comparable to Accelerating Rate Calorimeter tests in a fraction of the time it would take to run a Heat-Wait-Search test.

- Scanning Mode
  - Constant temperature range
  - Constant power
  - Fire exposure
- Isothermal
- Adiabatic (Heat-Wait-Search)

### Key Features of the MMC 274 *Nexus*<sup>TM</sup>

- Internal sample heater
- Stirring of the sample
- Injection of liquid sample
- Venting of the sample bomb
- Pressure measurement



Remote Calorimeter Options:

- Humidity-controlled container
- Refrigerated container
- Explosion-proof container

# MMC 274 Nexus™ – Application Examples

## MMC 274 Nexus™ operating like a DSC.




Epoxy adhesive at constant power input to sample



Pressure signal gives more information on solvents and decomposition compared to conventional DSC.

## MMC 274 Nexus™ operating like a an Accelerating Rate Calorimeter.

5 g Vazo67 solution in toluene. 10 ml titanium sample bomb fully compensated, to  $\Phi=1.0$  using VariPhi™ heater.



# Technical Data

The MMC 274 Nexus™ can be used for a wide range of applications in different areas, e.g.:

- Chemical Processing Industries
  - Process development and safety
  - Kinetics, □Hr
  - Vapor pressures, VLE
- Storage and Transportation of Chemicals
- Physical Properties Measurements
- Batteries
  - Thermal stability of whole cells
  - Component testing (cathode materials, electrolytes)
- Energetic Materials Testing

Parameter	MMC Operating Mode		
	Scanning	H-W-S	Isothermal
Temperature Scanning Rate (K/min)	up to 2	up to 2	NA
Sample Size (g)	0.1 to 5	0.1 to 5	0.1 to 5
Sample Volume (ml)	0.1 to 6	0.1 to 6	0.1 to 6
Sensitivity Heat Flow (μW/g)	250	200	25
Sensitivity (temp. rate) (K/min)	0.02		
Dynamic Range (W)	0.00006 to 0.5	0.0001 to 350	0.0000125 to ND
Temperature Range (°C)	25 to 500		
Pressure Range (bar)	0 to 150		
Sample Composition	Multi-Phase		
Adiabatic Tracking Rate - Fast module (K/min)	up to 200		
Adiabatic Tracking Rate - Standard module (K/min)	up to 25		
Options	Injection, Stirring, High Temp. Application Specific Containers		

**Exchangeable calorimeter module**



Technical specifications subject to change

Leading Thermal Analysis .

**NETZSCH**

**NETZSCH-Gerätebau GmbH**

Wittelsbacherstraße 42 · 95100 Selb, Germany  
 Phone: +49 9287 881-0 · Fax: +49 9287 881-505  
 e-mail: at@netzsch.com

www.netzsch.com